

MRI COMPATIBILITY EXABLATE NEURO


INSIGHTTEC

The **Exablate Neuro** is compatible with the following GE and Siemens MR models:

GE® MRI Scanners*

(Minimum SW requirements)

- SIGNA Artist 1.5T (DV26)
- Optima MR450w/GEM 1.5T (DV24)
- SIGNA HDxt 1.5T (HD23)
- SIGNA Premier 3T (Rx28)
- SIGNA Architect 3T (DV26)
- Discovery MR750 & 750w/GEM 3T (DV24)
- SIGNA HDxt 3T (HD15)



Siemens Healthineer MRI Scanners*

(Minimum SW requirements)

- MAGNETOM Altea 1.5T (XA20)
- MAGNETOM Sola 1.5T (XA20)
- MAGNETOM Aera 1.5T (VE11E)
- MAGNETOM Vida 3T (XA20)
- MAGNETOM Lumina 3T (XA20)
- MAGNETOM Prisma 3T (VE11E)
- MAGNETOM Prisma Fit 3T (VE11E)
- MAGNETOM Skyra 3T (VE11E)

*All Siemens and some GE scanners require the Exablate Compatibility Package from MR vendor for installation. Not all MR software versions are compatible with Exablate Neuro. Please confirm with your local Insightec representative.

